

RANDOX

1 SAMPLE • 1 TEST • 10 INFECTIONS

POINT OF CARE TESTING

Sexually Transmitted Infection Array

The Vivalytic enables sample to answer, cartridge based molecular diagnostics at the point of care. Powered by Randox Biochip Technology, Vivalytic provides the broadest cartridge based STI testing in the market. Detecting 10 bacterial, viral and protozoan infections, the STI test provides a comprehensive sexual health profile.

For more information contact: marketing@randox.com



 **BOSCH** 

RANDOX

SEXUALLY TRANSMITTED INFECTIONS (STI)

Sample Type: Swab or Urine

Sample Volume: 300 μ L

Detection Method: Randox Biochip Technology (end-point PCR)

Time to Result: 2 hours 30 minutes

INFECTIONS DETECTED

Chlamydia trachomatis (CT)

Neisseria gonorrhoeae (NG)

Trichomonas vaginalis (TV)

Mycoplasma genitalium (MG)

Treponema pallidum (Syphilis) (TP)

Herpes simplex virus 1 (HSV-1)

Herpes simplex virus 2 (HSV-2)

Haemophilus ducreyi (HD)

Mycoplasma hominis (MH)

Ureaplasma urealyticum (UU)

For more information contact: marketing@randox.com

