

# PathTROL™ Allergy Control Sera

## Multi-Analyte Allergy Controls

PathTROL™ Allergy Controls are liquid control sera, containing allergen specific IgE at two levels that are designed for the continuous, long term internal quality control of allergy tests in clinical and research laboratories.

PathTROL™ Allergy Controls contain specific human IgE to:

- d01 House Dust Mite (*D. pteronyssinus*)
- d02 House Dust Mite (*D. farinae*)
- e01 Cat Epithelia
- f17 Hazelnut
- f13 Peanut
- f02 Milk
- g06 Timothy Grass (*P. pratense*)
- g12 Cultivated Rye (*S. cereale*)
- w06 Mugwort (*A. vulgaris*)

The analytes are present in their native form in a physiological serum matrix, which are comparable to patients serum.

### Eurobio-Scientific

Eclipse House  
Curtis Road  
Dorking  
Surrey

Phone: 01306 888777  
Fax: 01306 883883  
E-mail: [sales@eurobioscientific.com](mailto:sales@eurobioscientific.com)  
Web: [www.eurobio-scientific.co.uk](http://www.eurobio-scientific.co.uk)

- Independent controls, suitable for most allergy assays and instruments
- Ready to use, liquid controls
- PathQAS™ Allergy peer-to-peer iQC scheme

PW80101	PATHTROL ALLERGY CONTROL SERA - LEVEL 1	5x1ml
PW80201	PATHTROL ALLERGY CONTROL SERA - LEVEL 2	5x1ml