

# PathTROL™ Tryptase Positive Controls

## High & Low Level Positive Controls

PathTROL™ Tryptase Controls are lyophilised control plasmas, containing Human  $\alpha$ -Tryptase and  $\beta$ -Tryptase at two levels that are designed for the continuous long term internal quality control of Tryptase tests in clinical and research laboratories.

- Independent controls, suitable for use with most Tryptase assays
- At two levels: Positive (20-30ng/ml) and High Positive (80-100ng/ml)
- PathQAS™ Tryptase, peer-to-peer iQC scheme included

### Eurobio Scientific

Eclipse House  
Curtis Road  
Dorking  
Surrey

Phone: +44 (0)1306 888777  
Fax: +44 (0)1306 883883

Web: [www.eurobio-scientific.com](http://www.eurobio-scientific.com)

PW85105 PathTROL™ Tryptase Positive Control 5x0.5ml  
PW85205 PathTROL™ Tryptase High Positive Control \* 5x0.5ml

*Also available; PathTROL™ Multi-Analyte Allergy Control sera*

PW80101 PathTROL™ Allergy Control - Level 1 5x1ml  
PW80201 PathTROL™ Allergy Control - Level 2 5x1ml